



diamedica

the market leader in anaesthetic equipment for challenging environments

Vaporiser

The Diamedica Vaporiser is a low resistance vaporiser suitable for both drawover and continuous flow anaesthesia. Thermally stable, the vaporiser gives a consistent output over a wide temperature range.

This is a versatile low resistance vaporiser:

- suitable for drawover or continuous flow anaesthesia
- available in Halothane/Isoflurane, or Sevoflurane versions

It requires no maintenance unless thymol from Halothane builds up. This can be cleared within seconds using fresh Halothane. It is used successfully and reliably in over 50 countries worldwide.

- Stainless steel construction
- Scale 0 to 5%
- Capacity 150ml
- Weight 2.6Kg dry
- Keyfill option available

Training and Support

Diamedica's commitment to our customers is:

- To provide, free of charge to all customers, ongoing after sales support via email or telephone for as long as the customer has their equipment
- To provide training for medical and technical personnel to enable them to maximise the potential of the equipment safely and economically

This vaporiser is supplied as part of the Diamedica Glostavent®, ECO2, and Portable anaesthesia machines.

Size: 16 x 10 x 15 cm

Weight: 2.6 kg



Diamedica (UK) Limited, Grange Hill Industrial Estate,
Bratton Fleming, Barnstaple, Devon, EX31 4UH, UK

Tel: +44 (0)1598 710066 Fax: +44 (0)1598 710055
Email: info@diamedica.co.uk www.diamedica.co.uk

© Copyright Diamedica (UK) Ltd